


Filter set 49 488049 - 9901

excitation: G 365
beamsplitter: FT 395
emission: BP 445/50


DAPI



Filter set 47 HE 489047 - 9901

excitation: BP 436/25
beamsplitter: FT 455
emission: BP 480/40


CFP



Filter set 38 HE 489038 - 9901

excitation: BP 470 / 40
beamsplitter: FT 495
emission: BP 525 / 50


eGFP



Filter set 46 HE 489046 - 9901

excitation: BP 500/25
beamsplitter: FT 515
emission: BP 535/30

YFP




Filter set 43 HE

489043 - 9901

excitation: BP 550 / 25
beamsplitter: FT 570
emission: BP 605 / 70

TRITC / Cy3.




Filter set 50

488050 - 9901

excitation: BP 640/30
beamsplitter: FT 660
emission: BP 690/50

cy5.



Filter set 63 HE

489063 - 0000

excitation: BP 572/25
beamsplitter: FT 590
emission: BP 629/62

mRFP



(typical curves)

